

LEG & ARM LIFT

The preparation of a surgical site on either the arms or legs of a patient often requires lifting, holding and lowering procedures. This job is normally undertaken by operating room staff, and can be physically demanding.

This unique electrical lift uses a hand control to lift the patient's leg or arm into the correct position. The limb can also be lowered using the hand control to ensure a smooth and safe process.

The lift operates by battery, and can easily be re-charged over-night.

The leg lift can accommodate a patient in a supine, abdominal or sitting position.

Article no.: 2009-01

Leg lift, electrical

Article no.: 2006-90

Lifting arm for 2 legs or abdominal pos., 90cm



Opitek International

Tel.: +45 39661644 • Fax: +45 39661655 • mail@opitek-international.dk

www.opitek-international.dk