

All glasses are made with the highest quality SCHOTT glass and offer .75mm lead equivalency. Most styles are available in Single, Bifocal, and Progressive prescriptions.



ITEM: LT100 **Ultralite**

Manufactured of lightweight nylon and designed to provide protection for the entire eye area. The curved front and temple bar with rubber tips allow for a close secure fit. Available in Black, Silver, Red, and Carbon Gray. Weight: 58 grams



ITEM: LT500 **Maxi**

The Maxi is slightly smaller than the Ultralite and includes a soft rubber nose bridge for increased comfort. Available in Black, Silver, Red, and Navy. Weight: 63 grams



ITEM: LT200 **Euro**

A standard safety frame which offers excellent protection and large viewing area. Universal form fit bridge. Weight: 75 grams



ITEM: LT300 **Astro**

These flexible safety frames can be adjusted for a more custom fit. The front of the frame is fitted snugly to the brow to help alleviate fluid run-off into the eyes. Weight: 75 grams



ITEM: LT5050 **Wrap**

A lightweight, wrap-around nylon frame with built-in adjustable nose piece to distribute weight evenly over the nose. The earpieces are attached by spring hinges designed to withstand the rigors of hospital use. Weight: 85 grams



ITEM: LT600 **MetalFlex**

These frames are constructed of optical steel base metal electroplate. Each pair is made with an adjustable non-slip silicone nose pad for maximum comfort. Weight: 62 grams



ITEM: LT400 **Fitover**

The Fitover is designed for wear over prescription glasses. The oversized frame offers excellent protection and large viewing area. Weight: 92 grams

