



Lightweight. Proven protection.™

PRODUCT CATALOG



RADIATION PROTECTION PRODUCTS

Quality you can trust

Lite Tech, Inc. was founded in 1992, with the acquisition of the patented XENOLITE business and technology from DuPont. Our products are developed from the longstanding DuPont technology and reputation of lightweight, eco-friendly, radiation protection innovation.

Headquartered near Philadelphia, Lite Tech is the only fully integrated x-ray protection manufacturer with an in-house R&D Lab, engineers, and ANSI representation on the IEC committee. Lite Tech produces its own proprietary core/sheet materials, which gives greater control over quality and standards compliance. XENOLITE products are backed by a thorough, vigorous testing process which meets the stringent world standard IEC 61331-1 (both 80 kV & 100 kV).

First in producing lightweight material, first in manufacturing a lead free option, and first with bi-layer technology, Lite Tech has and continues to develop unrivalled leading-edge capability in the combined sciences of radiation physics and complex highly-filled plastics processing. XENOLITE was the world's first lightweight, multi-element, composite x-ray protection material designed to replace lead-vinyl while offering a significant weight reduction.

Lite Tech products are sold in over 24 countries and at every interaction with each of our clients, we demonstrate our commitment to delivering the most innovative products to the highest standards of quality and compliance.

TO ORDER CALL (800) 647-5483
OUTSIDE THE U.S. CALL (610) 650-8690

Core Materials

Lead Free

XENOLITE® 800-NL — a lead-free, super-lightweight 2-element composite. XENOLITE NL is recyclable and safe for non-hazardous disposal.

PRODUCT CODE NUMBER IS 800

XENOLITE® 700-TB — a lead-free, lightweight 2-element composite. XENOLITE TB is recyclable and safe for non-hazardous disposal.

PRODUCT CODE NUMBER IS 700

XENOLITE® 400-S (Strata) — a lead-free, lightweight material featuring bi-layer radiation protection technology. XENOLITE Strata is recyclable and safe for non-hazardous disposal.

PRODUCT CODE NUMBER IS 400

Light-Weight Lead Composite

XENOLITE® 600-B — the original lightweight composite, 3-element material, with some lead content. XENOLITE is mid-priced and recyclable.

PRODUCT CODE NUMBER IS 600

Lead

EVAL 900-Pb — a 100% lead-based material. EVAL is used for lower-priced, non-weight sensitive products as well as gamma-radiation (Nuclear Medicine) applications. This material is also used for x-ray procedures (CT-assist, thorax) over 120 kV.

PRODUCT CODE NUMBER IS 900

Colors and Customization

Lite Tech offers a wide range of colors and fabrics for its XENOLITE aprons. To see the full selection of fabric options, visit www.xenolitexray.com. All orders are custom-made and can be embroidered with your logo or monogram.



Vest / Skirt Aprons

All XENOLITE Vest / Skirt Aprons create maximum weight distribution between the shoulders and hips which eliminates stress on the upper and lower back. Vest / Skirt Aprons also provide greater flexibility to the wearer with regard to sitting, bending, or stooping. The skirt is designed for complete overlap to provide maximum protection.

Vest / Skirt sizes can be mixed to provide maximum comfort and fit.

All Vest/Skirt Aprons come with standard features:

- Shoulder pads
- Pockets on vest and skirt
- Hanging loops on skirt


To determine sizing, please see page 15.



ITEM: 620VS / 720VS / 820VS

Weight Relief Vest / Skirt Apron

- Single-piece, front vest panel makes this model extremely light
- Velcro closure at right side and right shoulder of vest
- Unisex sizes



ITEM: 630VS / 730VS / 830VS

Vest / Skirt Apron

- The original design in wrap around protection
- Overlapping panels on vest allow for greater protection as well as greater adjustment
- Male and female sizes



ITEM: 625VS / 725VS / 825VS

Full Overlap Vest / Skirt Apron

- Full overlapping front panels allow for weight reduction and maximum protection
- Front panel closes at right side and right shoulder with Velcro
- Male and female sizes

ITEM: 635VS / 735VS / 835VS

Elastic Back Saver Vest / Skirt Apron

- Elastic Velcro straps cross at lower back to provide excellent lumbar support
- Eliminates weight on shoulders
- Velcro closure at left shoulder
- Male and female sizes



Full Wrap Aprons

All XENOLITE Full Wrap Aprons offer maximum protection with superior comfort and are designed to help reduce upper back stress during long procedures.

All Wrap Aprons come with standard features:

- Shoulder pads
- Pocket
- Adjustable 7" (18 cm) elastic belt comes standard with the Special Procedure Apron (6201 / 7201 / 8201) and the Overwrap Apron (6250 / 7250 / 8250)

To determine sizing, please see page 15.



ITEM: 6200TS / 7200TS / 8200TS

Tabbard Style Apron

- Right shoulder and side Velcro closure allow for easy access
- Simple, single-piece design makes this model extremely light
- Unisex sizes



ITEM: 6201 / 7201 / 8201

Special Procedure Apron

- Overlapping panels for greater protection and adjustment
- Adjustable 7" (18 cm) elastic belt distributes weight on hips and shoulders
- Velcro side seams allow for separation when bending or sitting without compromising protection
- Male and female sizes



ITEM: 6240 / 7240 / 8240

Diamond Back Apron

- One-piece, wrap-around opening in back with attached Lumbar support
- Front overlap eliminated for lighter weight
- Male & Female Sizes



ITEM: 6250 / 7250 / 8250

Overwrap Apron

- Full, overlapping, front panels allow for weight reduction and maximum protection
- Front panel closes at right side and right shoulder with Velcro
- Adjustable 7" (18 cm) elastic belt distributes weight on hips and shoulders
- Male and female sizes

ITEM: 6320 / 7320 / 8320

Quick Release Full Wrap Apron

Worn over the scrub suit and under the O.R. gown for quick removal without breaking the sterile field after x-ray procedures are completed

- Unisex Sizes



Frontal Aprons

All XENOLITE Frontal Aprons allow for freedom of movement, maximum flexibility, and optimal comfort.

The open back design makes all these aprons exceptionally light.

All Frontal Aprons come with standard features:

- Shoulder pads
- Pocket
- Male/Female sizes

To determine sizing, please see page 15.



ITEM: 610E / 710E / 810E

Elastic Tab Apron

- Elastic tabs allow for comfortable fit around the waist
- Reduces shoulder stress
- Body-hugging design shifts weight from shoulders to waist



ITEM: 610A / 710A / 810A

Adjustable Weight Apron

- Crisscross shoulder design keeps apron in place when bending or stooping
- Velcro tabs make it easy to put on and take off
- Reduces bunching on back and on shoulders
- Optional Velcro strap closure and/or 6" built-in elastic back panels



ITEM: 610K / 710K / 810K

Back Support Apron

- Unique design distributes weight between shoulders and lower body
- Adjustable 7" (18 cm) elastic support belt offers greater lower back support
- Belt attaches to side with key ring
- Velcro shoulder closure



ITEM: 610R / 710R / 810R

Shoulder Release Apron

- Crisscross, elastic back provides greater weight distribution
- Velcro tabs at both shoulders and waist allow for quick removal
- 4" Velcro strips for better adjustment

Frontal Aprons

All XENOLITE Frontal Aprons allow for freedom of movement, maximum flexibility, and optimal comfort.

The open back design makes all these aprons exceptionally light.

All Frontal Aprons come with standard features:

- Shoulder pads
- Pocket
- Male/Female sizes

To determine sizing, please see page 15.



ITEM: 610NV / 710NV / 810NV

No Velcro Apron

- 3" (7.5cm) wide elastic tabs allow for comfortable fit around the waist with buckle closure
- Body-hugging design shifts weight from shoulders to waist



ITEM: 610B / 710B / 810B

Buckle Apron

- Features a 2" (5 cm) wide adjustable waist strap with buckle
- Velcro shoulder closure
- Optional Velcro strap closure and/or 6" built in elastic back panels





ITEM: 610S / 710S / 810S

Strap Apron

- Features a cotton strap closure which is double stitched for durability
- Velcro shoulder closure

ITEM: 610LP / 710LP / 810LP

Lateral Protection Apron

All of the comfort advantages of our frontal aprons, plus additional coverage when turning to the side



ITEM: 610Q / 710Q / 810Q

O.R. Quick Drop Apron

Worn over the scrub suit and under the O.R. gown for quick removal without breaking the sterile field after x-ray procedures are completed



ITEM: 640T / 740T / 840T

Thyroid Collar

Available in Standard (fits most) or Large size (over 20"/51cm neck). Additional styles available.



ITEM: 650H / 750H / 850H

Humeral Shield



ITEM: LT15403

LiteTouch Sterile Attenuation Gloves

These lead- and powder-free thin gloves provide a super sensitive touch.



ITEM: 640H / 740H / 840H

Half Aprons

Velcro closure system; 0.50 mm Pb equivalency

Item	U.S. Size	Metric Size
X40H-1	8x10"	20x25
X40H-2	12x14"	30x35
X40H-3	16x18"	40x45
X40H-4	18x24"	45x60
X40H-5	Set of 3 8x10"	Set of 3 20x25
	12x14"	30x35
	16x18"	40x45

ITEM: 5200

Carry Bag

This durable bag is constructed of 500 Denier Nylon and carries up to two aprons. Color: Black

Monogramming is available for this item.



ITEM: 5100

Elastic Back Support Belt

Adjustable 7" (18 cm) elastic belt connects to apron with key ring.






ITEM: 9880

Mobile Apron Rack

Using only 24" of floor space, this mobile rack is perfect for rooms where floor space is limited. The 8 prongs can hold up to 12 aprons and thyroid collars.

Specifications:

48" tall. Color: White



ITEM: 9890

Valet Apron Rack

Constructed of sturdy painted metal and features wheels that lock in place, this valet holds up to 300 lbs of aprons. Five non-removable hangers are included with each unit. Additional hangers are available for purchase.

Specifications:

54" high, 36" wide, and 24" deep. Color: White



ITEM: 5136

8 Peg Apron Rack

Constructed of steel with rubber tips, this rack holds four aprons and thyroid collars. There are four mounting holes.

Specifications:

3" high, 24" wide. Color: White



ITEM: 9870

Z Frame Apron Rack

Our 4x4, heavy-duty, mobile, nesting Z-rack offers superior stability for storage of lead xray aprons. Large, heavy-duty, locking, swivel casters allow the rack to be easily moved. Five non-removable hangers are included with each unit.

Specifications:

55" high x 51" wide x 23" deep. Base Colors: white, black, & silver.

Protective Eyewear

All glasses are made with the highest quality SCHOTT glass and offer .75mm lead equivalency. Most styles are available in Single, Bifocal, and Progressive prescriptions.



ITEM: LT100 **Ultralite**

Manufactured of lightweight nylon and designed to provide protection for the entire eye area. The curved front and temple bar with rubber tips allow for a close secure fit. Available in Black, Silver, Red, and Carbon Gray. Weight: 58 grams



ITEM: LT500 **Maxi**

The Maxi is slightly smaller than the Ultralite and includes a soft rubber nose bridge for increased comfort. Available in Black, Silver, Red, and Navy. Weight: 63 grams



ITEM: LT200 **Euro**

A standard safety frame which offers excellent protection and large viewing area. Universal form fit bridge. Weight: 75 grams



ITEM: LT300 **Astro**

These flexible safety frames can be adjusted for a more custom fit. The front of the frame is fitted snugly to the brow to help alleviate fluid run-off into the eyes. Weight: 75 grams



ITEM: LT5050 **Wrap**

A lightweight, wrap-around nylon frame with built-in adjustable nose piece to distribute weight evenly over the nose. The earpieces are attached by spring hinges designed to withstand the rigors of hospital use. Weight: 85 grams



ITEM: LT600 **MetalFlex**

These frames are constructed of optical steel base metal electroplate. Each pair is made with an adjustable non-slip silicone nose pad for maximum comfort. Weight: 62 grams



ITEM: LT400 **Fitover**

The Fitover is designed for wear over prescription glasses. The oversized frame offers excellent protection and large viewing area. Weight: 92 grams



U.S. Sizing

VEST SIZE	WxL	SUIT SIZE	DRESS SIZE	CHEST	MAX. HEIGHT
X-Small	20x20"	34-35	2-4	32/33"	5'8"
Small	21x21"	36-37	6-8	34/35"	5'8"
Medium	22x22"	38-40	10-12	36/39"	5'10"
Large	23x23"	42-44	14-16	40/42"	6'0"
X-Large	24x24"	46-48	18-20	44/46"	6'2"
XX-Large	25x25"	50-52	22-24	48/50"	6'3"

SKIRT SIZE	WxL	SUIT SIZE	DRESS SIZE	CHEST	MAX. HIPS
X-Small	24x19"	34-35	2	21/24"	30"
Small	28x20"	36-37	4-10	25/28"	34"
Medium	32x21"	38-40	10-14	29/32"	38"
Large	36x22"	42-44	16-18	33/36"	42"
X-Large	40x23"	46-48	20-22	37/40"	46"
XX-Large	44x24"	50-52	24-26	41/44"	50"

WRAP SIZE	WxL	SUIT SIZE	DRESS SIZE	CHEST	MAX. HEIGHT
X-Small	20x34"	34-35	2-4	32/33"	5'6"
Small	21x36"	36-37	6-8	34/35"	5'8"
Medium	22x38"	38-40	10-12	36/39"	6'0"
Tall-Medium	22x40"	38-40T	10-12	36/39"	6'0"
Large	23x40"	42-44	14-16	40/42"	6'0"
Tall-Large	23x43"	42-44T	14-16	40/42"	6'3"
X-Large	24x40"	46-48	18-20	44/46"	6'0"
XX-Large	25x42"	50-52	22-24	48/50"	6'2"

MENS FRONTAL APRON SIZES	WxL	SUIT SIZE	CHEST	MAX. HEIGHT
Small	24x36"	36-37	34/35"	5'8"
Medium	24x38"	38-40	36/39"	6'0"
Large	24x40"	42-44	40/42"	6'0"
Tall-Large	24x43"	42-44T	40/42"	6'3"
X-Large	25x40"	46-47	44/45"	6'0"
XX-Large	26x42	48-50	46/48"	6'2"

WOMENS FRONTAL APRON SIZES	WxL	SUIT SIZE	CHEST	MAX. HEIGHT
X-Small	20x34"	2-4	32/33"	5'6"
Small	22x34"	4-6	34/35"	5'6"
Medium	22x36"	8-10	36/39"	5'8"
Tall-Medium	22x40"	8-10	36/39"	6'0"
Large	24x40"	14-16	40/42"	6'0"
X-Large	25x40"	18-20	43/46"	6'1"

Metric Sizing

VEST SIZE	W (cm)	SUIT SIZE	DRESS SIZE	CHEST
LENGTHS 53 / 55 / 60 / 63				
X-Small	51	44	34-36	80-84 cm
Small	53	46	36-38	85-90 cm
Medium	56	48-50	38-40	91-98 cm
Large	58	52-54	42-46	99-107 cm
X-Large	61	56	48-50	108-115 cm
XX-Large	64	58	52	116-130 cm

NOTE: Specify male or female

SKIRT SIZE	W (cm)	SUIT SIZE	DRESS SIZE	WAIST	MAX. HIPS
LENGTHS 50 / 55 / 60 / 63					
X-Small	61	44	34-36	53-62 cm	77 cm
Small	71	46	36-38	63-72 cm	87 cm
Medium	81	48-50	38-40	73-82 cm	97 cm
Large	91	52-54	42-44	83-92 cm	107 cm
X-Large	101	56	46-48	93-102 cm	117 cm
XX-Large	112	58	50-52	103-112 cm	127 cm

NOTE: Specify male or female

WRAP SIZE	W (cm)	SUIT SIZE	DRESS SIZE	CHEST	MAX. HIPS
LENGTHS 95 / 100 / 105 / 110 / 115 / 120					
X-Small	51	44	34-36	80-84 cm	77 cm
Small	53	46	36-38	85-90 cm	87 cm
Medium	55	48-50	38-40	91-96 cm	97 cm
Large	58	52-54	42-44	99-107 cm	107 cm
X-Large	61	56	46-48	108-115 cm	117 cm
XX-Large	63	58-60	50-52	116-130 cm	127 cm

NOTE: Specify male or female


MENS FRONTAL APRON SIZES	W (cm)	EURO SUIT SIZE	UK SUIT SIZE	CHEST
LENGTHS 100 / 105 / 110 / 115 / 120				
Small	61	46	36	86-90 cm
Medium	61	48-50	38-40	91-98 cm
Large	61	52-54	42-44	99-107 cm
X-Large	63	56-58	44-48	108-115 cm
XX-Large	66	58-62	48-54	116-130 cm

WOMENS FRONTAL APRON SIZES	W (cm)	EURO SUIT SIZE	UK SUIT SIZE	CHEST
LENGTHS 95 / 100 / 105 / 110 / 115				
X-Small	51	34-36	6-8	80-84 cm
Small	56	36-38	10-12	85-89 cm
Medium	56	40-42	14-16	90-94 cm
Large	61	42-44	18-20	95-107 cm
X-Large	63	46-50	20-22	108-115 cm

Determining size

Before placing an order, take body measurements. Use a cloth measuring tape to get accurate measurements.

1. Measure around chest at underarm
2. Measure around waist
3. Measure around hip at widest point
4. Measure length from top of shoulder to length required



To learn more and view a sizing video, visit

xenolitexray.com



975 Madison Avenue
Norristown, PA 19403
T: (610) 650-8690 / (800) 647-5483
F: (610) 650-8694
info@xenolitexray.com
www.xenolitexray.com

CE
0120