

MEDICAL EQUIPMENT FOR MRI

HARLOFF
650 Ford Street
Colorado Springs, CO 80915

Harloff's MR-Conditional Carts are built primarily of non-ferrous aluminum, plastic, and stainless steel for safe daily use in highly magnetic environments.

Now Harloff offers two cabinet sizes – Standard Size and Narrow Body.

With packages for Anesthesia and Emergency requirements, multiple drawer configurations, and extra MR-Conditional accessories, Harloff is the right choice for storage, security and delivery needs in the Magnetic Resonance Imaging Suite.



MRI LABELING INFORMATION

The following products were determined to be MR-Conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005:

Cart 1, Model: MR7K-MAN, 7 Drawer Keylocking MR Cart with Anesthesia Accessory Package

Cart 2, Model: MR6B-EMG, 6 Drawer Breakaway Locking MR Cart with Emergency Cart Package

Cart 3, Model: MRN7K, 7 Drawer Narrow Body Keylocking MR Cart (tested with Emergency Accessory Package, #EMG)

This product included the following:

- MR-4TILT – 4 Tilt Bin Accessory with Mounting Clip
- MR-5TILT – 5 Tilt Bin Accessory with Mounting Clip
- MR-6TILT – 6 Tilt Bin Accessory with Mounting Clip
- MR-3GSHARPS – 3 Gallon Sharps Container


Non-clinical testing demonstrated that each product is MR-Conditional and can be used in the MRI environment according to the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient magnetic field of 720-Gauss/cm or less

IMPORTANT NOTE: These products are intended for use inside of the MRI environment (e.g., in the MR system room, close to the scanner). However, they should not be utilized directly inside of the MR system (e.g., inside of the bore of the scanner).



**MR-CONDITIONAL CARTS
AT 3-TESLA OR LESS**



**MEDICAL EQUIPMENT
FOR THE MRI ENVIRONMENT**

HARLOFF
www.harloff.com

Contract # V797P-4938A
05/15/2006 through 05/14/2011
GSA Advantage!

MAGNETIC RESONANCE EMERGENCY CARTS

MR-Conditional at 3-Tesla or Less

STANDARD SIZE BODY



MR6B

Six Drawer, Emergency Cart
Standard Cabinet, Standard Package
45.75" H x 30.5" W x 23.75" D
(116.2 cm H x 77.5 cm W x 60.3 cm D)
Net Weight: 126 lbs (57.2 kg)

MR6B-EMG

Six Drawer, Emergency Cart
Standard Cabinet, Specialty Package
45.75" H x 34.5" W x 25.5" D
(116.2 cm H x 87.6 cm W x 64.8 cm D)
With IV Pole: 68"/174.7 cm H
Net Weight: 144 lbs (65.4 kg)

MR7B

Seven Drawer, Emergency Cart
Standard Cabinet, Specialty Package
45.75" H x 30.5" W x 23.75" D
(116.2 cm H x 77.5 cm W x 60.3 cm D)
Net Weight: 126 lbs (57.2 kg)

MR7B-EMG

Seven Drawer, Emergency Cart
Standard Cabinet, Standard Package
45.75" H x 34.5" W x 25.5" D
(116.2 cm H x 87.6 cm W x 64.8 cm D)
With IV Pole: 68"/174.7 cm H
Net Weight: 144 lbs (65.4 kg)

STANDARD SIZE BODY



MR6K

Six Drawer, Anesthesia Cart
Standard Cabinet, Standard Package
45.75" H x 30.5" W x 23.75" D
(116.2 cm H x 77.5 cm W x 60.3 cm D)
Net Weight: 126 lbs (57.2 kg)

MR6K-MAN

Six Drawer, Anesthesia Cart
Standard Cabinet, Specialty Package
45.75" H x 32" W x 25" D
(116.2 cm H x 81.3 cm W x 63.5 cm D)
With Back Rail: 68.5"/174 cm H
Net Weight: 139 lbs (63.1 kg)

MR7K

Seven Drawer, Anesthesia Cart
Standard Cabinet, Standard Package
45.75" H x 30.5" W x 23.75" D
(116.2 cm H x 77.5 cm W x 60.3 cm D)
Net Weight: 126 lbs (57.2 kg)

MR7K-MAN

Seven Drawer, Anesthesia Cart
Standard Cabinet, Specialty Package
45.75" H x 32" W x 25" D
(116.2 cm H x 81.3 cm W x 63.5 cm D)
With Back Rail: 68.5"/174 cm H
Net Weight: 139 lbs (63.1 kg)

NARROW BODY



MRN6B

Six Drawer, Narrow Body
Emergency Cart, Standard Package
46" H x 20" W x 23.5" D
(116.8 cm H x 50.8 cm W x 59.7 cm D)
Net Weight: 85 lbs (38.6 kg)

MRN6B-EMG

Six Drawer, Narrow Body
Emergency Cart, Specialty Package
46" H x 20" W x 23.5" D
(116.8 cm H x 50.8 cm W x 59.7 cm D)
With IV Pole: 68.75"/174.6 cm H
Net Weight: 113 lbs (51.4 kg)

MRN7B

Seven Drawer, Narrow Body
Emergency Cart, Standard Package
46" H x 20" W x 23.5" D
(116.8 cm H x 50.8 cm W x 59.7 cm D)
Net Weight: 85 lbs (38.6 kg)

MRN7B-EMG

Seven Drawer, Narrow Body
Emergency Cart, Specialty Package
46" H x 20" W x 23.5" D
(116.8 cm H x 50.8 cm W x 59.7 cm D)
With IV Pole: 68.75"/174.6 cm H
Net Weight: 113 lbs (51.4 kg)

NARROW BODY



MRN6K

Six Drawer, Narrow Body
Anesthesia Cart, Standard Package
46" H x 20" W x 23.5" D
(116.8 cm H x 50.8 cm W x 59.7 cm D)
Net Weight: 85 lbs (38.6 kg)

MRN6K-MAN

Six Drawer, Narrow Body
Anesthesia Cart, Specialty Package
68.5" H x 20" W x 23.5" D
(174 cm H x 50.8 cm W x 59.7 cm D)
Net Weight: 108 lbs (50.5 kg)

MRN7K

Seven Drawer, Narrow Body
Anesthesia Cart, Standard Package
46" H x 20" W x 23.5" D
(116.8 cm H x 50.8 cm W x 59.7 cm D)
Net Weight: 85 lbs (38.6 kg)

MRN7K-MAN

Seven Drawer, Narrow Body
Anesthesia Cart, Specialty Package
68.5" H x 20" W x 23.5" D
(174 cm H x 50.8 cm W x 59.7 cm D)
Net Weight: 108 lbs (50.5 kg)