



KANMED^o BABYBED

THE MODERN INTENSIVE CARE BABY BED THAT PROVIDES EFFICIENT ERGONOMICS.

KANMED^o BABYBED provides the best care for premature and newborn babies in medium and low intensive care situations. The cosy Baby Bed is developed together with neonatal nurses to guarantee a high level of ergonomics. The height of the Baby Bed can easily be electrically adjusted and the see-through side walls are easy to fold for the best access to the baby. The bed comes with a

nice designed tent to reduce the baby's exposure to light, air condition and sound. Most Kanmed Baby Beds are equipped with the KANMED^o BABYWARMER. Together they form a modern care station in Neonatal Intensive Care, Maternity and Delivery wards. KANMED^o BABYBED and KANMED^o BABYWARMER are found all over the world.

KANMED^o BABYBED

KANMED^o BABYBED IS AVAILABLE IN 3 VERSIONS:
(THE BABY BEDS ON THE PICTURES BELOW ARE EQUIPPED WITH THE KANMED^o BABYWARMER)



BB-1

KANMED^o BABYBED

Complete unit Includes:

Electrical height adjustment, Mattress, Tent with Tent Pole, Shelf, Power Cable, 19 cm high Foldable Side Walls (24 cm height is also available).



BB-3

KANMED^o BABYBED FIXED HEIGHT

Complete unit Includes:

Mattress, Tent with Tent Pole, Shelf, 19 cm high Foldable Side Walls (24 cm height is also available).



BB-4

KANMED^o BABYBED TWINS

Complete unit Includes:

Electrical height adjustment, Mattress, Tent with Tent Pole, Shelf, Power Cable, 24 cm high Foldable Side Walls (19 cm is also available). Made for 2 complete Kanmed Baby Warmers or one Baby Warmer with the twin water mattress and nest.

KANMED^o BABYWARMER

Do you know that the KANMED^o BABYWARMER is used in more than 1000 hospitals worldwide? KANMED^o BABYWARMER makes it possible for the premature baby to move from the incubator to a more baby friendly and growth oriented environment. It has been shown that the baby becomes calmer, sleeps better and grows faster on a warm water mattress in comparison with an incubator. The very important bonding between the parents and baby can start and the general care of the baby is easier in an open system. Also, the total cost of the Baby Warmer care is only a fraction of what incubator care costs. KANMED^o BABYWARMER is most often used in combination with the, CPAP and other and other vital signs monitoring.



KANMED^o BABYBED

OPTIONAL ACCESSORIES TO KANMED^o BABYBED



BB-10I DRAWER SYSTEM

Drawer system that is mounted on the back of the bed. Complete with 2 plastic store boxes, 13 litres each and a mounting bracket.



BB-XX BACK SHELF COMPLETE

Vertical poles 25mm, Horizontal rails 10x25 mm, 1 top shelf 500 x300 mm. Including brackets and screws.



BB-XX

Tents in different colours are available. Check colours and availability on www.kanmed.se



BB-XX BOTTLE HOLDER

Stainless steel bracket for holding hand disinfection bottles. Max dia 70-80.



BB-57 TENT KIT FOR COTS

Fits into a standard cot. Min width of cot xxcm. Can rotate to move tent away from cot. Is hold in place by the bed.



BB-XX CORNER BRACKET

Universal bracket for mounting for instance a vertical pole or other items.



BB-XX POWER CABLE

Exordium Narratio Humbugus

See www.kanmed.se for the latest information about products and accessories.

TECHNICAL DATA: KANMED BABYBED

ART. NO.	BB-1	BB-4	BB-3
DESCRIPTION:	Standard bed, Electrically adjustable height	Twin bed, Electrically adjustable height	Fixed height bed
INCLUDES:	Mattress, Tent, Tent pole, Shelf, 19 cm side walls (24 cm for twin version) power cable, manuals and export packaging.		
OUTER size:	760 mm x 560 mm	850 mm x 700 mm	760mm x 560mm
BED SURFACE size:	665 mm x 445 mm	765 mm x 640 mm	665mm x 445mm
Maximum load Bed surface:	10 Kg		
WEIGHT:	About 33 kg	About 40 kg	About 25 kg
WHEEL BASE:	700mm x 500mm		
WHEELS:	Diameter 125mm. Front wheels with brakes.		
TILT:	± 22°		
MATTRESS:	Polyurethane with PVC cover. (665x445x30mm)	Polyurethane with PVC cover. (765x640x30mm)	Polyurethane with PVC cover. (665x445x30mm)
TENT:	100% Cotton. Washable 60°C (all models)		
SHELF:	Min 400mm x 400mm, All models. Max load 10 Kg		
COLOUR:	Unisex light blue (all models)		
VOLTAGE:	100-240Vac, 50-60Hz	100-240Vac, 50-60Hz	NA
POWER CONSUMPTION:	127 VA	127VA	NA
APPLIED STANDARDS: All beds carries the CE-Mark	SS EN 60 60I. Conforms to MDD	SS EN 60 60I. Conforms to MDD	
PROTECTION:	Class I	Class I	Class I

MAIN ACCESSORIES:	
BB-20 FOAM MATTRESS, BB-48 Foam mattress for twin bed	
BB-21/I TENT YELLOW/YELLOW (Standard tent)	BB-22 Twin bed tent
BB-102 BACK SHELF COMPLETE	Universal corner bracket
BB-101 DRAWER SYSTEM	Bottle holder
700-0721 Front wall, 700-0729 Front wall 24 cm height	700-0732 Front wall twin
700-0723 Side wall, 700-0731 Side wall 24 cm height	700-0728 Side wall twin
700-0722 Rear panel, 700-0733 Rear panel 24 cm height	700-0733 Rear panel twin
700-0750 Y-cable	
BB-103 TENT+TENT POLE HOLDING DEVICE PACK OF 5 PC.	