

SINGLE-USE, STERILE TOURNIQUET CUFFS

Single-use, Sterile Tourniquet Cuffs (complete with stocking sleeve padding) are the latest addition to our range, providing sterility when the risk of infection is key.

Because they are sterile, these cuffs can be used in front of sterile drapes, giving the surgeon optimum operating space – particularly useful in paediatrics or for procedures such as total knee replacement or ulna nerve decompression.








If a towel slips or is breached, these cuffs do not compromise the sterile field, giving added security against cross infection.



Single-use, Sterile Tourniquet Cuffs
providing sterility when the risk of infection is key



Single-use Sterile Tourniquet Cuffs (box of 10)

Description	Code
 Single-use Tourniquet Cuff – Sterile – Paediatric	40800
 Single-use Tourniquet Cuff – Sterile – Adult Arm	40820
 Single-use Tourniquet Cuff – Sterile – Small Leg	40840
 Single-use Tourniquet Cuff – Sterile – Adult Leg	40860
 Single-use Tourniquet Cuff – Sterile – Conical Adult Leg	40868
 Single-use Tourniquet Cuff – Sterile – XL Adult Leg	40870
 Single-use Tourniquet Cuff – Sterile – XL Conical Adult Leg	40878

STERILE, SINGLE-USE

less time on cleaning / checking for re-use

SLIM, LOW PROFILE

designed for patient comfort

SECURE FASTENING

eliminates roll or twist

LIGHT WEIGHT

also latex free

FULL RANGE OF SIZES

colour coded for ease of recognition



DESIGN & INNOVATION



MADE IN THE UK



GLOBAL DISTRIBUTION