


Bonney Uterine Polyp Forceps

Product Specification

The Design Checklist

- ✓ Bonney Uterine Polyp Forceps have 6mm slender head sizes that improve visibility and access in comparison to the Rampley Sponge Holder, which has a 10mm
- ✓ Designed to the highest quality, used by the best clinicians
- ✓ Experience true control and clarity
- ✓ Designed to meet specified gynaecological requirements
- ✓ Slender strong jaw with precision aligned tips, suitable for many procedures such as polyp and coil removal
- ✓ Effective ratchet that enables a strong and robust grip, superior to the plastic Rampley Sponge Holder alternatives



Code	Pack size	Description
DT- BPF2001	20 units	Bonney Uterine Polyp Forceps 240mm, double packed