

Dental Syringe

Product Specification

The Design Checklist

- Suitable for use where needle stick injury is a concern (compliant with EU directive 2010/32/EU – prevention from sharps injuries)
- Popular top loading style to accommodate 2.2ml vials
- Barrel diameter is wide enough to take most vials
- Available in extra-long length for use with larger, high BMI patients
- Designed for administering anaesthetic during a wide range of procedures



Code	Pack size	Description
DT-DGS2001	20 units	Dental Syringe 85mm for 2.2ml vial, single packed
DT-DGSL2001	20 units	Dental Syringe Long 115mm for 2.2ml vial, single packed