

Episiotomy Scissors

Product Specification

The Design Checklist

- High quality cut
- It has a slender sharp double-bevelled blade for precise cutting
- Shorted flat blunt blade for support
- Sharp blades facilitate insertion in to the vagina



5.25mm tip width -
53mm tip length



67mm tip length

Code	Pack size	Description
DT-BSS2001	20 units	Braun Stadler Episiotomy Scissors Angular 145mm, double packed
DT-WPS2018	20 units	Waldmann Perineum Episiotomy Scissors, double packed