

Clearway Suction Handle

Product Specification

The Design Checklist

- Removes difficult foreign objects from ear and nose
- Ideal for A&E/ER and MIU
- Less traumatic for patient and family
- Can remove the need for ENT surgeon intervention or theatre time
- UK manufactured
- 100% quality inspected
- Our innovative tip, has a soft malleable suctioning cap that forms an air-tight seal to securely grasp foreign objects
- 142° angled handle for improved visibility
- Suction control, particularly useful for smooth, rounded objects
- Precision manufactured to ensure consistent quality
- Universal suction attachment
- Readily available Bubble Tubing or Yankauer Tubing to be cut down at a suitable point once the handle has been removed



Code	Pack size	Description
DT-SVT1001-S	10 units	Clearway Suction Handle, single packed