

## Crocodile Fine Jaw Micro Forceps Product Specification

### The Design Checklist

- The best Single-use option available, making them the product of choice when removing Foreign Bodies from the Ear Cavity
- Commonly used for Wax Removal, Foreign Body Removal and Grommet Placement
- Precise 0.9 x 5mm jaw action that can be used with the smallest Specula
- For Paediatric use or where good vision is critical
- Smooth, matt black finish prevents any reflection when used with a Microscope providing you with enhanced clarity during procedures
- Optimal first-time use and accurate grip
- 100% inspection and dimensional checks on every Instrument



Code	Pack size	
DT-CFBF2001-S	20 units	Crocodile Fine Jaw Micro Forceps black finish, single packed