

**Fish Hook Retractor**

**Product Specification**

**The Design Checklist**

- Stainless steel hook, available in two sizes
- Sharp point for accurate positioning
- Latex free 2mm silicone 180mm ties
- Pack of five for varying wound sizes
- Fish Hook Retractor: 35mm Hook Length & 28mm Hook Width
- Paediatric Fish Hook Retractor: 28mm Hook Length & 22mm Hook Width



Code	Pack size	Description
DT-FHR0510	10 units	Fish Hook Retractors x 5 with Silicone Loop, double packed
DT-FHRP0510	10 units	Paediatric Fish Hook Retractors x 5 with Silicone Loop, double packed