


Keele Mucosectomy Forceps

Product Specification

The Design Checklist

- Ideal for use during Trans Oral Mucosectomy procedure
- Forceps are used to hold, control, manipulate and stabilise tissue in the Oral Cavity and manoeuvre the tongue during Cauterisation
- Used for removing Carcinomas at the back of the Tongue
- 140° angled tip for enhanced access
- 1×2 Teeth to maintain firm specimen retention



Code	Pack size	Description
DT-KMF1025	10 units	Keele Mucosectomy Forceps, 250mm, 16mm tip length, 140° angle