

### Oral Biopsy Forceps

### Product Specification

#### The Design Checklist

- Elongated shaft to aid procedure
- Two cup sizes allowing you flexibility to take a better sample
- Single-use Oral Biopsy Forceps eliminates cross contamination, inadequate samples & repeat procedures
- Ensures a precise biopsy is taken to assist diagnosing oral cancer or any other abnormality in the oral and laryngeal mucosa



**DTR Medical®**

AN  
INNOVIA  
MEDICAL  
COMPANY

| Code         | Pack size | Description   |
|--------------|-----------|---|
| DT-OBF1001-S | 10 units  | Oral Biopsy Forceps 90mm stainless steel, single packed |
| DT-OBF1001-S | 10 units  | Oral Biopsy Forceps Large 120mm stainless steel, single |