

## Thudichum Nasal Speculum

## Product Specification

### The Design Checklist

- Thudichum Nasal Speculum featuring a strong flexible stainless steel construction with smooth edges
- Balanced sprung blades for ease of use
- Multipack containing three individual sizes
- Particularly suitable for clinic use
- New and improved size 0 for paediatric use



Code	Pack size	Description
DT-TNS2003	10 x 3 units	Thudichum Nasal Specula, multipack (sizes 1, 2 & 3), single packed
DT-TNS2000	20 units	Thudichum Nasal Specula size 0 (Extra Small), single packed
DT-TNS2011	20 units	Thudichum Nasal Specula size 1 (Small), single packed
DT-TNS2022	20 units	Thudichum Nasal Specula size 2 (Medium), single packed
DT-TNS2033	20 units	Thudichum Nasal Specula size 3 (Large), single packed