



Double Ended Punctum Dilator Product Specification

The Design Checklist

- Manufactured in the UK to ensure the instrument provides consistent precision every time
- Used in Punctoplasty to dilate the Lacrimal Punctum for syringing and flushing
- Designed with high quality polished tips to avoid trauma during entry to the punctum
- Knurled handle supports improved grip which allows increased accuracy in dilation procedures
- Well balanced and double ended to suit your procedural needs
- 17mm and 21mm tapered ends for increased range of dilation



Ø0.3mm tip



Ø0.4mm tip



Code	Pack size	Description
DT-PDL2014	20 units	Double Ended Punctum Dilator 140mm