

Barron Suction Handles

Product Specification

The Design Checklist

- Complete range of sizes from 3 to 7Fg with suction control
- Working length of 77mm



Code	Pack size	Description
DT-BSH2003	20 units	Barron Suction Handle 1mm (3Fg), double packed
DT- BSH2004	20 units	Barron Suction Handle 1.35mm (4Fg), double packed
DT- BSH2005	20 units	Barron Suction Handle 1.67mm (5Fg), double packed
DT- BSH2006	20 units	Barron Suction Handle 2mm (6Fg), double packed
DT- BSH2007	20 units	Barron Suction Handle 2.3mm (7Fg), double packed