

**Capsulorhexis Forceps**

**Product Specification**

**The Design Checklist**

- Curved shaft offers greater flexibility
- Utrata style tips offer improved functionality



CXF2001- 0.95mm tip length x 0.4mm tip width



CXFF2003- 0.6mm tip length x 0.3mm tip width



Code	Pack size	Description
DT-CXF2001	20 units	Capsulorhexis Forceps, Utrata Tips, 12mm Curved Shafts, Round Handles
DT-CXFF2003	20 units	Capsulorhexis Fine Forceps, 0.3mm Tips