

Frazier Suction Handles

Product Specification

The Design Checklist

- Strong construction suitable for bending to shape during procedure
- Single-use avoids the need to re-straighten or risk breakage
- Tip options include non-fenestrated, fenestrated with two holes and cross fenestrated with four holes
- Expanded range with 150° degree bend
- NEW Pre-bent handles
- NEW 5mm (15Fg) variants now available




With suction control



Cross fenestrated



Without suction control



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

FRAZIER SUCTION HANDLES

2mm (6Fg)

FSH2006	Frazier Suction Handle 200 x 2mm OD (6Fg) straight, double packed	20 units
FSH2006-B	Frazier Suction Handle 190 x 2mm OD (6Fg) 150° bend, double packed	20 units
FSF2006	Frazier Suction Handle 200 x 2mm OD (6Fg) straight, fenestrated, double packed	20 units
FSF2046	Frazier Suction Handle 200 x 2mm OD (6Fg) straight, cross fenestrated, double packed	20 units
FSN2426	Frazier Suction Handle 184 x 2mm OD (6Fg), 42° bend, without suction control, double packed	20 units

2.66mm (8Fg)

FSH2008	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, double packed	20 units
FSH2008-B	Frazier Suction Handle 190 x 2.66mm OD (8Fg) 150° bend, double packed	20 units
FSF2008	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, fenestrated, double packed	20 units
FSF2048	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, cross fenestrated, double packed	20 units

3mm (9Fg)

FSH2009	Frazier Suction Handle 200 x 3mm OD (9Fg) straight, double packed	20 units
FSH2009-B	Frazier Suction Handle 190 x 3mm OD (9Fg) 150° bend, double packed	20 units
FSF2009	Frazier Suction Handle 200 x 3mm OD (9Fg) straight, fenestrated, double packed	20 units
FSF2049	Frazier Suction Handle 200 x 3mm OD (9Fg) straight, cross fenestrated, double packed	20 units
FSN2429	Frazier Suction Handle 184 x 3mm OD (9Fg), 42° bend, without suction control, double packed	20 units

4mm (12Fg)

FSH2012	Frazier Suction Handle 200 x 4mm OD (12Fg), 150° bend, double packed	20 units
FSF2012	Frazier Suction Handle 200 x 4mm OD (12Fg), fenestrated, 150° bend, double packed	20 units
FSF2412	Frazier Suction Handle 200 x 4mm OD (12Fg), cross fenestrated, 150° bend, double packed	20 units
FSN24212	Frazier Suction Handle 184 x 4mm OD (12Fg), 42° bend, without suction control, double packed	20 units