


## Kerrison Rongeur

## Product Specification

### The Design Checklist

- NEW 100mm variant with 90 degrees up bite complements our existing DCR Range that includes our Punctum Dilator, Lacrimal Probes, Hartmann Nasal Speculum, Kilner “Catspaw” Retractor, and Suction Handles.
- 200mm shaft length with 40 degrees up bite
- Well-balanced stainless-steel construction
- Specifically designed for cutting and removing small pieces of bone during surgery
- Eliminates the risk of cross contamination between procedures



Code	Pack	Description
DT-KRU1101	10 units	Kerrison Rongeur Punch 1mm jaw, length 200mm with 40 degrees up bite, double packed
DT-KRU1201	10 units	Kerrison Rongeur Punch 2mm jaw, length 200mm with 40 degrees up bite, double packed
DT-KRU1301	10 units	Kerrison Rongeur Punch 3mm jaw, length 200mm with 40 degrees up bite, double packed
DT-KRU1401	10 units	Kerrison Rongeur Punch 4mm jaw, length 200mm with 40 degrees up bite, double packed
DT-KRU1501	10 units	Kerrison Rongeur Punch 5mm jaw, length 200mm with 40 degrees up bite, double packed
DT-KRUF11029	10 units	Kerrison Rongeur Punch 2mm wide, length 100mm with 90 degrees up bite, double packed
DT-KRUF11039	10 units	Kerrison Rongeur Punch 3mm wide, length 100mm with 90 degrees up bite, double packed