

**Aspirating Needle**

**Product Specification**

**The Design Checklist**

- Aspirating Needle with Stainless steel cannula with 18g needle for laparoscopic aspiration
- Single-use avoids difficult to clean cannula and needle with stepped joint
- Extra long 335mm length suitable for patients of all sizes
- Optimum first-time sharpness of needle for ease of use
- Available in 3.5 and 5mm diameter cannulas complete with Luer connection



Code	Pack size	Description
DT- ASN2018	20 units	Aspirating Needle 335 x 3.5mm (18g) with Luer, double packed
DT- ASN5018	50 units	Aspirating Needle 335 x 3.5mm (18g) with Luer, double packed