

## Adson Forceps


## Product Specification

### The Design Checklist


- Narrow tips to effectively grip and hold delicate tissue in place
- Flexibility across a wide range of procedures



Adson Toothed Forceps  
- 2mm tip width



Adson Non-Toothed Forceps - 15mm tip length x 2mm tip width



Adson Forceps with Plateau- 1.4mm tip width



Fine Tip Adson Toothed Forceps - 22mm tip length x 1.5mm tip width



**DTR Medical**®

AN INNOVIA  
MEDICAL  
COMPANY

Code	Pack size	Description
DT-ADS5001-S	50 units	Adson Non-Toothed Forceps 120mm, single packed
DT-ATS5001-S	50 units	Adson Toothed Forceps 120mm, single packed
DT-ATP2012-D	20 units	Adson Forceps 1x2 with plateau 120mm, double packed
DT-ATP2015-D	20 units	Adson Forceps 1x2 with plateau 150mm, double packed
DT-ATS2012F	20 units	Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm, double packed