

GENERAL & VASCULAR SURGERY



DTR Medical[®]

AN
INNOVIA
MEDICAL
COMPANY

Sterile Single-use
Enhancing surgical workflow



About Innovia Medical

Innovia Medical is a family of specialty surgical companies that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com

Contents

Forceps <i>Crile Artery, Mosquito, Kocher Artery, Kocher-Ochsner, Adson, McIndoe, Fine Suture, Debakey, Gillies, Baby Overholt, Allis Tissue, Rampley Sponge Holders, Spencer Wells & Strassburg Forceps</i>	4-7
Needle Holders <i>Mayo-Hegar, Crilewood, Halsey, Ryder, Kilner & DeBakey Needle Holders</i>	8-9
Scissors <i>Kilner, Iris, Vannas, Mayo, Reynolds, Metzenbaum, Stevens Tenotomy & SuperCut Scissors</i>	9-14
Clamps <i>Backhaus Towel & Tubing Clamps</i>	15
Bipolar Forceps <i>Standard, Adson, Cawthorne, Bipolar Cable, Monopolar Forceps & Bayonet Bipolar Forceps</i>	16-18
Curettes <i>Beckmann Adenoid Curette</i>	19
Clamp Covers <i>Sterile & Non-Sterile Clamp Covers</i>	19
Silicone Slings	20
General & Vascular <i>Heparin Flushing Needles</i>	20
Tibbs Arterial Cannula	21
Phlebectomy Hooks	22
Dissectors <i>Watson Cheyne Dissector</i>	22
Retractors <i>Kilner Alae Hook Retractor</i>	22
Suction Handles <i>Bellucci, Long, Micro Laryngeal, Bronchoscopy, Oesophageal, Sigmoidoscope, Frazier, Micro, Suction Regulator, Poole & Yankauer</i>	23-26
Tip Shields	26
Tubing Products <i>Silicone & Sterile Silicone Tube Sets</i>	27-30
Others <i>Silicone Bands, Blue Clamps & Scalpel Blade Handles</i>	31
Scope Warmer	32
Instrument Wipes	32

Forceps

CRILE ARTERY FORCEPS

CFS2001-D	Crile Artery Forceps 145mm Straight	20 units
CFC2001-D	Crile Artery Forceps 145mm Curved on Flat	20 units

- Commonly used for clamping blood vessels or tissue
- Effective locking ratchet and robust grip



Straight



Curved



Crile Artery Forceps Curved


MOSQUITO FORCEPS

MQT5001-S	Mosquito Forceps 125mm Straight	50 units
MQC2012	Micro Mosquito Forceps 120mm Curved	20 units
MQC2025-D	Mosquito Forceps Curved 125mm	20 units
MQC2125-D	Mosquito Forceps Curved Toothed 1x2 125mm	20 units

- Frequently used for small blood vessels
- Type of Hemostatic Forceps commonly used to control bleeding
- Robust locking ratchet and serrated jaws for holding tissue in place



Straight



Curved



Micro Mosquito Curved 120mm


KOCHER ARTERY FORCEPS

KAF2014-D	Kocher Artery Forceps 140mm	20 units
-----------	-----------------------------	----------

- Effective locking ratchet and robust grip
- Type of Hemostatic Forceps commonly used to control Bleeding
- 1x2 teeth at the tip, our Forceps ensure a secure grasp on heavy Tissue and when clamping off bleeding Blood Vessels



37mm tip length




Kocher Artery Forceps 140mm

Forceps

KOCHER-OCHSNER FORCEPS

KOF2181 Kocher-Ochsner 1x2 Toothed Forceps Straight 185mm 20 units


- Secure grasp on heavy tissue and when clamping off bleeding blood vessels



ADSON FORCEPS

ADS5001-S	Adson Non-Toothed Forceps 120mm	50 units
ATS5001-S	Adson Toothed Forceps 120mm	50 units
ATP2012-D	Adson Forceps 1x2 with plateau 120mm	20 units
ATP2015-D	Adson Forceps 1x2 with plateau 150mm	20 units
ATS2012F	Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm	20 units


- Narrow tips to effectively grip and hold delicate tissue in place
- Flexibility across a wide range of procedures




Adson Toothed Forceps



Adson Non-Toothed Forceps



Adson Forceps with Plateau




Fine Tip Adson Toothed Forceps

McINDOE FORCEPS

MDF5001-D McIndoe Dissecting Forceps 150mm, Non-Toothed 50 units

- Fine, non-toothed serrated jaws are ideal when re-positioning delicate Tissue
- Can be used to hold wound dressings and packing




Forceps


FINE SUTURE FORCEPS

SPF2001-S	Fine Suture Forceps Plain 105mm	20 units
SSF2001-S	Fine Suture Forceps Serrated 105mm	20 units

- Provides flexibility required in delicate procedures
- Deliver an excellent grip to help secure tissue and suture materials



Plain



Serrated



Fine Suture Forceps Serrated

DEBAKEY FORCEPS

DAF2016-D	DeBakey Atraumatic Dissecting Forceps 160mm	20 units
-----------	---	----------


- Used to grasp soft tissue and organs during surgery
- Specially designed teeth maintain the grip on the tissues, but do not tear



GILLIES FORCEPS

GDF2001-D	Gillies Toothed Dissecting Forceps 160mm	20 units
-----------	--	----------


- Available with fine, narrow jaws
- 1:2 serrated teeth
- Less pressure required to adequately grip harder material such as skin



BABY OVERHOLT FORCEPS

BOC2014	Baby Overholt Artery Forceps Curved, 140mm	20 units
---------	--	----------

- Curved to allow easy access to veins




Forceps

ALLIS TISSUE FORCEPS

ATF2341 Allis Tissue Forceps 3:4 teeth, 160mm 20 units

- Enables a strong grasp in surgical procedures
- Finger ring, ratchet forceps
- Ensure an effective hold of slippery tissue



RAMPLEY SPONGE HOLDERS

RSH2181 Rampley Sponge Holder 180mm 20 units

RSH2241 Rampley Sponge Holder 240mm 20 units

- Strong jaws with an effective locking ratchet for effectively holding swabs and gripping tissue during procedures
- Robust grip



SPENCER WELLS FORCEPS

SWF2010 Spencer Wells Artery Forceps 120mm Straight 20 units

SWF2015 Spencer Wells Artery Forceps 150mm Straight 20 units

SWF2018 Spencer Wells Artery Forceps 180mm Straight 20 units

SWF2020 Spencer Wells Artery Forceps 200mm Straight 20 units

- Used during surgery to compress the Artery, clamp, and seal small Blood Vessels or hold the Artery out of the way
- Gripping jaws and finely serrated blades



Needle Holders

MAYO-HEGAR NEEDLE HOLDERS

MNH2013	Mayo-Hegar Needle Holder 130mm	20 units
MNH2015	Mayo-Hegar Needle Holder 150mm	20 units
MNHG2013-D	Mayo-Hegar Needle Holder with groove 130mm	20 units

- Engineered with serrated jaws that allow for an effective grip
- Central groove to ensure the Needle can be securely held



MNH2015



MNHG2013-D




Mayo, 150 mm

Mayo with groove, 130 mm

NEEDLE HOLDERS

CNH2015	Crilewood Needle Holder 150mm	20 units
HNH2013	Halsey Needle Holder 130mm	20 units
RNH2015-D	Ryder Needle Holder 150mm	20 units

- Ensure a firm grip when Suturing
- Robust locking ratchet




CNH2015



HNH2013



RNH2015-D



Crilewood Needle Holder

Halsey Needle Holder


Ryder Needle Holder

Needle Holders

NEEDLE HOLDERS

KNH2013	Kilner Needle Holder 135mm	20 units
DNH2018-D	DeBakey Needle Holder 180mm	20 units

- Designed with narrow, atraumatic teeth to ensure a secure grasp on fine sutures and delicate tissue whilst minimising damage



KNH2013



DNH2018-D



Kilner Needle Holder 135mm





DeBakey Needle Holder 180mm

Scissors

KILNER SCISSORS

KSS2001-S	Kilner Scissors 115mm	20 units
-----------	-----------------------	----------

- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips
- Ensure ease when cutting and dissecting delicate tissue in surgical procedures




Scissors


IRIS SCISSORS

ISS2111-S	Iris Scissors Straight 90mm	20 units
ISC2111-S	Iris Scissors Curved 90mm	20 units
ISS2110-D	Iris Scissors Straight 110mm	20 units
ISC2110-D	Iris Scissors Curved 110mm	20 units

- Enhanced accuracy when it comes to delicate surgeries
- First time sharpness which cut even at the tips
- Narrows to a fine and sharp point



Straight



Curved




NEW

SUPERCUT IRIS SCISSORS

SISS2111	SuperCut Iris Scissors Straight 90mm	20 units
SISC2111	SuperCut Iris Scissors Curved 90mm	20 units


- One serrated and one sharp cutting blade allows a better grip and cut than traditional alternatives
- Sharp tips make the Scissors ideal for precise cutting of soft tissue, sutures, and wound dressings
- Available with straight or curved blades

IMAGE COMING SOON




Vannas Scissors


- High quality stainless steel cross action scissors with fine sharp tips
- Precise blade action
- Curved option offers greater flexibility for delicate procedures
- Ideal for your most delicate practice, such as the removal of fine Sutures
- Equipped with a ridged handle and a double leaf spring mechanism for self-opening action, ensuring enhanced control



Straight



Curved



VANNAS SCISSORS


VSS2081	Vannas Scissors Straight 85mm	20 units
VSS2121	Vannas Scissors Straight 125mm	20 units
VSC2081	Vannas Scissors Curved 85mm	20 units

Scissors


MAYO SCISSORS

MSS2141	Mayo Scissors Straight Bevelled Blade 140mm	20 units
MSC2141	Mayo Scissors Curved Bevelled Blade 140mm	20 units
MSS2171	Mayo Scissors Straight Bevelled Blade 170mm	20 units
MSC2171	Mayo Scissors Curved Bevelled Blade 170mm	20 units

- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips
- Choice of curved or straight bevelled blades



Straight



Curved



Mayo Scissors Straight 170mm





Mayo Scissors Curved 170mm

REYNOLDS SCISSORS



RSC2014-D	Reynolds Scissors Curved 145mm	20 units
-----------	--------------------------------	----------

- Fine tipped scissors enable a precise action when cutting and dissecting soft tissue
- 145mm length helps keep the clinician's hand out of view in demanding procedures



Metzenbaum Scissors

- First time performance for cutting delicate tissue and blunt dissection
- Specially engineered to be fine and precise



METZENBAUM SCISSORS

ZSC2121-D	Baby Metzenbaum Scissors Curved 120mm	20 units
ZSC2141-D	Metzenbaum Scissors Curved 145mm	20 units
ZSC2161-D	Metzenbaum Scissors Curved 160mm	20 units
ZSS2181-D	Metzenbaum Scissors Straight 180mm	20 units
ZSC2181-D	Metzenbaum Scissors Curved 180mm	20 units
ZSS2201-D	Metzenbaum Scissors Straight 200mm	20 units
ZSC2201-D	Metzenbaum Scissors Curved 200mm	20 units

NEW

SUPERCUT BABY METZENBAUM SCISSORS

SZCS2121	SuperCut Baby Metzenbaum Scissors 120mm, Curved	20 units
----------	---	----------

- Ideal for fine cutting of soft tissue and blunt dissection
- One serrated blade and one sharp cutting blade allows a superior grip and cut





IMAGE COMING SOON

Scissors

STEVENS TENOTOMY SCISSORS

TSS2001-S	Stevens Tenotomy Scissors 115mm Straight	20 units
TSC2001-S	Stevens Tenotomy Scissors 115mm Curved	20 units

- Stainless steel sharp/sharp Scissors
- First time sharpness which cut even at the tips
- Smooth and elegant long handle
- Allows increased grip and control when cutting in small or narrow areas



NEW

SUPERCUT STEVENS TENOTOMY SCISSORS

STSS2001	SuperCut Stevens Tenotomy Scissors 115mm Straight	20 units
STSC2001	SuperCut Stevens Tenotomy Scissors 115mm Curved	20 units

- Ideal for delicate cutting and dissecting soft tissue
- One serrated blade and one sharp cutting blade allows a superior grip and cut
- Available with straight or curved blades






IMAGE COMING SOON

Clamps

CLAMPS

BTC2013-D	Backhaus Towel Clamp 135mm	20 units
TCF2051	Tubing Clamps 150mm	20 units

- Ensures durability and are manufactured using the highest quality stainless steel
- Secure grip for clamping tubing without damage





i You may also be interested in Silicone Tubing see p27-30

Bipolar Forceps

STANDARD BIPOLAR FORCEPS

SU770513044	Standard Bipolar Forceps 130 x 0.5mm, 3m, 4mm Plugs	20 units
SU770515044	Standard Bipolar Forceps 150 x 0.5mm, 3m, 4mm Plugs	20 units
SU770520044	Standard Bipolar Forceps 200 x 0.5mm, 3m, 4mm Plugs	20 units
SU771015044	Standard Bipolar Forceps 150 x 1mm, 3m, 4mm Plugs	20 units
SU771520044	Standard Bipolar Forceps 200 x 1.5mm, 3m, 4mm Plugs	20 units
SU772015044	Standard Bipolar Forceps 150 x 2mm, 3m, 4mm Plugs	20 units
SU772020044	Standard Bipolar Forceps 200 x 2mm, 3m, 4mm Plugs	20 units

- Fine stainless steel tips with a range of sizes from 0.5mm - 2mm
- Adjustable and universal 4mm plug that is tested to comply with the latest 'Flying Leads' UK legislation
- High quality cables and connectors deliver consistently repeatable results



ADSON BIPOLAR FORCEPS


SU790515044	Adson Bipolar Forceps 150 x 0.5mm, 3m, 4mm Plugs	20 units
SU791015044	Adson Bipolar Forceps 150 x 1mm, 3m, 4mm Plugs	20 units



Bipolar Forceps

CAWTHORNE BIPOLAR FORCEPS

SU760508044	Cawthorne Bipolar Forceps 80 x 0.5mm, 3m, 4mm Plugs	20 units
SU761008044	Cawthorne Bipolar Forceps 80 x 1mm, 3m, 4mm Plugs	20 units
SU760511044	Cawthorne Bipolar Forceps 110 x 0.5mm, 3m, 4mm Plugs	20 units
SU761011044	Cawthorne Bipolar Forceps 110 x 1mm, 3m, 4mm Plugs	20 units




BIPOLAR CABLE

BPC5030	Bipolar Cable, 3.0m lead, 4mm pins UK/EU Fitting	50 units
---------	--	----------



MONOPOLAR FORCEPS


SU671220044	Monopolar Forceps 200mm Straight x 1.2mm tip serrated	20 units
-------------	---	----------



Bayonet Bipolar Forceps

- Fine stainless steel tips with 1mm or 2mm options
- Adjustable and universal 4mm plug that is tested to comply with the latest 'Flying Leads' UK legislation
- High quality cables and connectors deliver consistently repeatable results

SU781020044



SU782020044



BAYONET BIPOLAR FORCEPS

SU781020044	Bayonet Bipolar Forceps 200 x 1mm, 3m, 4mm Plugs	20 units
SU782020044	Bayonet Bipolar Forceps 200 x 2mm, 3m, 4mm Plugs	20 units

Curettes

BECKMANN ADENOID CURETTE

BAC2010	Beckmann Adenoid Curette, 10mm Curette Edge, 220mm	20 units
BAC2013	Beckmann Adenoid Curette, 13mm Curette Edge, 220mm	20 units
BAC2015	Beckmann Adenoid Curette, 15mm Curette Edge, 220mm	20 units
BAC2017	Beckmann Adenoid Curette, 17mm Curette Edge, 220mm	20 units

- Sharp and precise scoops to enable an effective Adenoidectomy
- Serrated handle to provide a secure grip whilst scraping Tissue



Clamp Covers


STERILE CLAMP COVERS

CCR2031	Clamp Cover Red (Small) 50mm	20 units
CCW2041	Clamp Cover White (Med) 120mm	20 units
CCB2051	Clamp Cover Blue (Large) 120mm	20 units

CLAMP COVERS (NON-STERILE)

SM3040	White - Medium	5 x 1m
SM3050	Blue - Large	5 x 1m

- Soft silicone jaw covers ensure good grip with reduced trauma
- Custom sterile packs cut to preferred lengths available on request



CCR2031



SM3040



SM3050

Silicone Slings

SILICONE SLINGS

1.5mm Slings

SM2100	Multicoloured Slings 1.5mm packed x 5	20 units
SM2103	Red Slings 1.5mm	25 units
SM2104	Blue Slings 1.5mm	25 units
SM2105	Yellow Slings 1.5mm	25 units
SM2106	White Slings 1.5mm	25 units

2mm Slings

SM2002	Multicoloured Slings 2mm packed x 5	20 units
SM2008	Sling Multipack 2mm, (Red x 3, Blue x 3, Yellow x 1, White x 1)	20 units
SM2030	Red Slings 2mm	25 units
SM2040	Blue Slings 2mm	25 units
SM2050	Yellow Slings 2mm	25 units
SM2060	White Slings 2mm	25 units

- Atraumatic oval cross section for less damage during surgery
- Radiopaque and biocompatible medical grade silicone
- Colour coded for your specific procedure:
 - Blue: Vein
 - Red: Artery
 - Yellow: Nerve
 - White: Tendon



General and Vascular

HEPARIN FLUSHING NEEDLES

HFN2001	Heparin Flushing Needle 48mm	20 units
HFN2081	Heparin Flushing Needle 85mm	20 units


- Available in two distinct lengths to enable better access to the Vein
- Single-use removes the risk of cross-contamination making it very suitable for your high-risk cases



UNIQUE RANGE

Tibbs Arterial Cannula

- Unique range of Single-Use Tibbs Arterial Cannula with a high-quality stainless-steel construction
- 8 variants are provided, available in Micro, Small, Large, X Large, and XX Large
- The Cone and Olive variants increase flexibility, ensuring you get the correct size for the specific vessel when flushing them with heparin solution
- UK manufactured




TIBBS STANDARD SIZES

TAC20LC	Tibbs Arterial Cannula Large Cone, 48mm	20 units
TAC20SC	Tibbs Arterial Cannula Small Cone, 48mm	20 units
TAC20LO	Tibbs Arterial Cannula Large Olive, 48mm	20 units
TAC20SO	Tibbs Arterial Cannula Small Olive, 48mm	20 units



TIBBS SPECIAL SIZES

TAC20MO	Tibbs Arterial Cannula Micro Olive 2.5mm OD, 48mm	20 units
TAC20XLC	Tibbs Arterial Cannula X Large Cone 8mm OD extra-long, 69mm	20 units
TAC20XXLC	Tibbs Arterial Cannula XX Large Cone 12.25mm OD extra-long, 69mm	20 units
TAC2003	Tibbs Arterial Cannula Set of 3 sizes (Large Olive, Large Cone & X Large Cone)	20 units



Phlebectomy Hooks

PHLEBECTOMY HOOKS

PBH2001	Phlebectomy Hook Size 1 (2mm)	20 units
PBH2002	Phlebectomy Hook Size 2 (1.5mm)	20 units
PBH2003	Phlebectomy Hook Size 3 (1mm)	20 units

- Allows for an effective removal of Varicose Veins that lie underneath the surface of the leg
- Handle shaped for comfort and control during procedures
- Readily available in 1.0, 1.5 and 2.0mm Hook sizes to ensure a precise incision





Dissectors

WATSON CHEYNE DISSECTOR

WCD2225	Watson Cheyne Dissector, 225mm	20 units
---------	--------------------------------	----------

- Flexible device can be used for the mobilisation and removal of degenerative Tissue in General Surgery
- Very fine blunt ends with slender curves which allow for high precision and control
- Avoids the use of traditional re-usable alternatives that are "knobby and fat"




Retractors

KILNER ALAE HOOK RETRACTOR

KHR2001	Kilner Alae Hook Retractor	20 units
---------	----------------------------	----------

- Effective retraction for surgical exposure on






Suction Handles

LONG SUCTION HANDLES

240mm		
BLSF2240	Micro Laryngeal Suction Handle 240 x 2mm OD (6Fg)	20 units
280mm		
BLS2280	Micro Laryngeal Suction Handle 280 x 2mm OD (6Fg)	20 units
330mm		
BLS2330	Micro Laryngeal Suction Handle 330 x 2mm OD (6Fg)	20 units
BLSF2339	Micro Laryngeal Suction Handle 330 x 3mm OD (9Fg)	20 units
360mm		
BLSF2369	Micro Laryngeal Suction Handle, 360 x 3mm OD (9Fg)	20 units
450mm		
BLS2450	Bronchoscopy Suction Handle 450 x 4mm OD (12Fg)	20 units
550mm		
BLS2550	Oesophageal Suction Handle 550 x 4mm OD (12Fg)	20 units
560mm		
SGH2001	Sigmoidoscope Suction Handle, 560 x 6mm OD (18Fg)	20 units

- Effective alternative to reprocessing very long narrow lumen
- Range of sizes from 240mm (6Fg) to 560mm (18Fg)





EXTENDED RANGE

Frazier Suction Handles

- Strong construction suitable for bending to shape during procedure
- Single-use avoids the need to restraighten or risk breakage
- Tip options include non-fenestrated, fenestrated with two holes and cross fenestrated with four holes
- Expanded range with 150° degree bend
- Pre-bent handles
- 5mm (15Fg) variant now available

Supplied straight* and pre bent
Bend to shape for your procedure

*Excluding 12Fg range



FRAZIER SUCTION HANDLES

2mm (6Fg)

FSH2006	Frazier Suction Handle 200 x 2mm OD (6Fg) straight	20 units
FSH2006-B	Frazier Suction Handle 190 x 2mm OD (6Fg) 150° bend	20 units
FSF2006	Frazier Suction Handle 200 x 2mm OD (6Fg) straight, fenestrated	20 units
FSF2046	Frazier Suction Handle 200 x 2mm OD (6Fg) straight, cross fenestrated	20 units
FSN2426	Frazier Suction Handle 184 x 2mm OD (6Fg), 42° bend, without suction control	20 units

2.66mm (8Fg)

FSH2008	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight	20 units
FSH2008-B	Frazier Suction Handle 190 x 2.66mm OD (8Fg) 150° bend	20 units
FSF2008	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, fenestrated	20 units
FSF2048	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, cross fenestrated	20 units

3mm (9Fg)

FSH2009	Frazier Suction Handle 200 x 3mm OD (9Fg) straight	20 units
FSH2009-B	Frazier Suction Handle 190 x 3mm OD (9Fg) 150° bend	20 units
FSF2009	Frazier Suction Handle 200 x 3mm OD (9Fg) straight, fenestrated	20 units
FSF2049	Frazier Suction Handle 200 x 3mm OD (9Fg) straight, cross fenestrated	20 units
FSN2429	Frazier Suction Handle 184 x 3mm OD (9Fg), 42° bend, without suction control	20 units

4mm (12Fg)

FSH2012	Frazier Suction Handle 200 x 4mm OD (12Fg), 150° bend	20 units
FSF2012	Frazier Suction Handle 200 x 4mm OD (12Fg), fenestrated, 150° bend	20 units
FSF2412	Frazier Suction Handle 200 x 4mm OD (12Fg), cross fenestrated, 150° bend	20 units
FSN24212	Frazier Suction Handle 184 x 4mm OD (12Fg), 42° bend, without suction control	20 units

5mm (15Fg)

FSH2015	Frazier Suction Handle 200 X 5mm OD (15Fg), 150° bend	20 units
---------	---	----------

Suction Handles

EXTENDED RANGE

MICRO SUCTION TUBES

86mm

MSV508113	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g)	50 units
MSV508113R	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g) with Suction Regulator	50 units
MSV508114	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g)	50 units
MSV508114R	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g) with Suction Regulator	50 units
MSV508116	Micro Suction (Verhoeven) Cream 86 x 1.65mm OD (16g)	50 units
MSV508117	Micro Suction (Verhoeven) White 86 x 1.47mm OD (17g)	50 units
MSV508119	Micro Suction (Verhoeven) Brown 86 x 1.05mm OD (19g)	50 units

110mm

MSV50114	Micro Suction (Verhoeven) Green 110 x 2.1mm OD (14g)	50 units
MSV50117	Micro Suction (Verhoeven) White 110 x 1.47mm OD (17g)	50 units
MSV50119	Micro Suction (Verhoeven) Brown 110 x 1.05mm OD (19g)	50 units
MSV50122	Micro Suction (Verhoeven) Black 110 x 0.7mm OD (22g)	50 units

120mm

MSV50214	Micro Suction (Verhoeven) Green 120 x 2.1mm OD (14g)	50 units
----------	--	----------

- Stainless steel tube with colour coded luer connector for microsuction use
- Also suitable for use with syringe as irrigation device
- Available in 86, 110 and 120mm lengths



UNIVERSAL ADAPTER

SUCTION REGULATOR

DSRS5001	Suction Regulator, Stainless Steel with Improved Grip	50 units
----------	---	----------

- House style regulator with control
- Luer connector suitable for Verhoeven and other micro suction tubes



Suction Handles

SPECIAL SIZES

POOLE SUCTION HANDLES

PSHN1001	Poole Suction Handle Neonatal	10 units
PSH1001	Poole Suction Handle	10 units
PSHL1001	Poole Suction Handle Long	10 units

- Stainless steel for balance and precision
- Supplied with detachable outer to spread suction
- Available in three lengths 190, 240, 330mm to suit patients from neonates to high BMI adults




Poole Suction Handle

CHOICE OF ROSES

YANKAUER SUCTION HANDLES

YST1000	Yankauer Suction Handle 305mm with rose	10 units
YST1000R	Yankauer Suction Handle 305mm with modified rose	10 units

- Detachable rose with three (YST1000) or six holed modified rose (YST1000R)
- Ideal for tonsillectomy



Yankauer Suction Handle 305mm with rose

Tip Shields

TIP SHIELDS

SM5020	Flat 13mm width Brown	100 units
SM5050	Round 6mm diameter Yellow	100 units
SM5060	Round 4mm diameter Red	100 units
SM5080	Round 2mm diameter Blue	100 units
SM5090	Round 1mm diameter Blue	100 units

- Supplied non-sterile
- Wide choice of colour, size, and shape
- Cost effective instrument protection
- Protects staff during reprocessing and sterilisation



Tubing Products

EASY TO USE SIZE GUIDE

SILICONE TUBING

SM135	ID 0.51	OD 0.94	15m roll
SM155	ID 0.64	OD 1.19	15m roll
SM205	ID 1.02	OD 2.16	15m roll
SM265	ID 1.57	OD 2.41	15m roll
SM285	ID 1.57	OD 3.18	15m roll
SM305	ID 1.98	OD 3.18	15m roll
SM335	ID 3.35	OD 4.65	15m roll
SM365	ID 3.18	OD 6.35	15m roll
SM405	ID 4.8	OD 7.9	15m roll
SM445	ID 6.4	OD 9.5	15m roll
SM465	ID 6.4	OD 11.1	15m roll
SM485	ID 6.4	OD 12.7	15m roll
SM525	ID 9.5	OD 12.7	15m roll
SM545	ID 9.5	OD 14.3	15m roll

- Supplied non-sterile
- Supplied in 15 metre packs for cutting to size
- Suitable for custom procedure trays



SUCTION HANDLE TUBING RECOMMENDATIONS

Tubing (non-sterile)	Tubing Diameters	Spigot	Sterile DSTS Code	Suitable for use Suction Handles
SM405	4.8mm ID x 7.9mm OD	4-6mm	DSTS40506	Aspirating Dissectors Bronchoscopy Frazier Magill Microlaryngeal Oesophageal Poole Sinus Suction Regulator
SM445	6.4mm ID x 9.5mm OD	7-11mm	DSTS44506	All of the above (exc Sinus) plus Barron, Schuknecht and Tapered

Sterile Single-Use Silicone Tube Sets

HOW TO USE: First 6 digits + order code followed by inner diameter (ID) and outer diameter (OD) then length (mm) e.g. DST205 + 1004 = DST2051004 (1.02 ID x 2.16 OD) 40mm

DSTS = Spigot	Inner Diameter	Outer Diameter	40	75	150	200	250	300	400	450	600	700	750-760	800	1000	1100	1350-1450	1500-1750	2000	2500	3600	
CODE	mm	mm																				
DST105	0.30	0.64						03														
DST205	1.02	2.16	10												10							
DST265									54													
DST285	1.57	3.18					1 or 2		03	04		07										
DST305	1.98	3.18				02																
DST325	2.65	4.88						03														
DST365	3.20	6.40												08		11				20		
DSTS365	3.20	6.40				02																
DST405	4.80	7.90			015			03			06		76-5		10		13	15 or 17				
DST415	5.00	8.30									06							15				
DSTS425	4.80	9.50											76									
DST445	6.40	9.50						03				07			10			15		25	36	
DSTS445	6.40	9.50									06				10							
DST465	6.40	11.10							04													36
DST545	9.50	12.70		007																		
DSTS545	9.50	14.30					02															
DST625	12.70	17.50																			25	
PACK SIZE KEY	10 (units)	20 (units)																				

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

Tubing Products

EXTENDED RANGE

STERILE SILICONE TUBE SETS

DST105

DST10503 Sterile Silicone Tube 0.3 ID x 0.64 OD x 300mm long 10 units

DST205

DST2051004 Sterile Silicone DCR Tube 1.02 ID x 2.16 OD x 40mm long 10 units

DST20510 Sterile Silicone Tube 1.02 ID x 2.16 OD x 1000mm long 10 units

DST265

DST26554 Sterile Silicone Tube 1.57 ID x 2.41 OD x 400mm long 10 units

DST285

DST2851 Sterile Silicone Tube 1.57 ID x 3.18 OD x 250mm long 10 units

DST2852 Sterile Silicone Tube 1.57 ID x 3.18 OD x 250mm long 2 pieces 10 units

DST28503 Sterile Silicone Tube 1.57 ID x 3.18 OD x 300mm long 10 units

DST28504 Sterile Silicone Tube 1.57 ID x 3.18 OD x 450mm long 10 units

DST28507 Sterile Silicone Tube 1.57 ID x 3.18 OD x 700mm long 10 units

DST305

DST30502 Sterile Silicone Tube 1.98 ID x 3.18 OD x 200mm long 10 units

DST325

DST32503 Sterile Silicone Tube 2.65 ID x 4.88 OD x 300mm long 10 units

DST365

DSTS36502 Sterile Silicone Tube 3.2 ID x 6.4 OD x 200mm long with Spigot 20 units


DST36508 Sterile Silicone Tube 3.2 ID x 6.4 OD x 800mm long 10 units

DST36511 Sterile Silicone Tube 3.2 ID x 6.4 OD x 1100mm long - House 20 units

DST36520 Sterile Silicone Tube 3.2 ID x 6.4 OD x 2000mm long 10 units

- Sterile for extra confidence where lengths over 270mm are impossible to validate
- Cost effective alternative to in-house production
- Full range of tube sizes cut to consistent lengths
- Tailored to suit needs including mixed packs
- Selected tube sets available with Spigot

Tube Dimensions



STERILE TUBING AND SPIGOT STARTER PACKS

Code	Description	Pack size
DSTS40506	Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 600mm long with Spigot	10 units
DSTS44506	Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 600mm long with Spigot	10 units

Tubing Products

STERILE SILICONE TUBE SETS (CONTINUED)

DST405

DSTS40506	Sterile Silicone Tube 4.8 ID x 7.9 OD x 600mm long with Spigot	10 units
DST405015	Sterile Silicone Tube 4.8 ID x 7.9 OD x 150mm long	10 units
DST40503	Sterile Silicone Tube 4.8 ID x 7.9 OD x 300mm long	10 units
DST40506	Sterile Silicone Tube 4.8 ID x 7.9 OD x 600mm long	10 units
DST40576-S	Sterile Silicone Tube 4.8 ID x 7.9 OD x 760mm long	10 units
DST40510	Sterile Silicone Tube 4.8 ID x 7.9 OD x 1000mm long	10 units
DST40513	Sterile Silicone Tube 4.8 ID x 7.9 OD x 1350mm long	20 units
DST40515	Sterile Silicone Tube 4.8 ID x 7.9 OD x 1500mm long	20 units
DST40517	Sterile Silicone Tube 4.8 ID x 7.9 OD x 1750mm long	20 units

DST415

DST41506	Sterile Silicone Tube 5.0 ID x 8.3 OD x 600mm long	10 units
DST41515	Sterile Silicone Tube 5.0 ID x 8.3 OD x 1500mm long	20 units

DST425

DSTS42576	Sterile Silicone Tube 4.8 ID x 9.5 OD x 760mm long with Spigot	10 units
-----------	--	----------

DST445

DSTS44506	Sterile Silicone Tube 6.4 ID x 9.5 OD x 600mm long with Spigot	10 units
DST44503	Sterile Silicone Tube 6.4 ID x 9.5 OD x 300mm long	10 units
DST44507	Sterile Silicone Tube 6.4 ID x 9.5 OD x 700mm long	10 units
DST44510	Sterile Silicone Tube 6.4 ID x 9.5 OD x 1000mm long	10 units
DSTS44510	Sterile Silicone Tube 6.4 ID x 9.5 OD x 1000mm long with Spigot	10 units
DST44515	Sterile Silicone Tube 6.4 ID x 9.5 OD x 1500mm long	20 units
DST44525	Sterile Silicone Tube 6.4 ID x 9.5 OD x 2500mm long	20 units
DST44536	Sterile Silicone Tube 6.4 ID x 9.5 OD x 3600mm long	20 units

DST465

DST46504	Sterile Silicone Tube 6.4 ID x 11.1 OD x 450mm long	20 units
DST46536	Sterile Silicone Tube 6.4 ID x 11.1 OD x 3600mm long	10 units

DST545

DST545007	Sterile Silicone Tube 9.5 ID x 14.3 OD x 75mm long	10 units
DSTS54502	Sterile Silicone Tube 9.5 ID x 14.3 OD x 250mm long with Spigot	10 units

DST625

DST62525	Sterile Silicone Tube 12.7 ID x 17.5 OD x 2500mm long	20 units
----------	---	----------

Others

SILICONE BANDS

DSB3211-S Silicone Bands 2.16 x 110mm diameter x 3 20 units

- Sterile for extra confidence
- Cost effective alternative to in-house production
- Latex free to replace elastic bands
- Customised options available



BLUE CLAMPS

DBC2051-S Blue Clamp Mediplast 20 x 1 unit
DBC5051-S Blue Clamp Mediplast 10 x 5 units



UNIVERSAL ADD ON

SCALPEL BLADE HANDLES

SBH2031 Scalpel Blade Handle No 3 20 units
SBH2041 Scalpel Blade Handle No 4 20 units

- Secure grip
- Eliminates cross infection



Scope Warmer

THEATRE EFFICIENCY

SCOPE WARMER

31SW10	ClearView Anti-Fog Scope Warmer, sterile packed for Single-use	10 units
31SW30	ClearView Anti-Fog Scope Warmer, sterile packed for Single-use	30 units

- Simple preheating to eliminate condensation
- Secures the surgeon's vision
- Warms full scope length to 37°C within 90 seconds, staying above body temperature for 120 minutes
- Long-acting to avoid disruption or repeated defogging



Instructional poster available
Contact us on +44 (0)1792 797910

Instrument Wipes

NETCELL® PVA INSTRUMENT WIPE

35-900	PVA Instrument Wipe	20 units
--------	---------------------	----------

- Quickly removes protein, blood and particles without damage to the instrument
- Manufactured from ultra-smooth, biocompatible PVA sponge that will not snag
- Soft, strong, flexible, fibre and snag free material



“Supporting visibility, quality and safety in your supply chain”**

Our instruments' packaging features the latest machine readable GS1-standard barcoded labels for product traceability.

The newly incorporated labelling technology is a timely development with an increasing number of health providers worldwide adopting an 'eProcurement' strategy and will provide machine readable product information on all labels to meet expected EU and other regulators demands.

Associated benefits include:

- Improve patient care and reduce costs
- Improve stock control
- One global language
- Accurately track equipment

“Unique device identification (UDI) compliance represents a significant development for our customers. Automation of information flow will have influences across the whole purchasing and supply process, providing the benefits of superior instrument traceability and improved patient safety.”

“We are working closely with Kodit UDI software to provide our labelling solutions to meet the latest international regulations and we serve as a reference site for other manufacturers to learn from.”

DTR Medical

For more information on our new product labelling, please contact Customer Services on +44 (0) 1792 797910 or dtrinfo@innoviamedical.com

** - GS1 – Global Traceability Standards

Matching reusable standards

We benchmark against reusable British and international standards:

BS 5194-2:1989

Specification for Instruments with Pivot Joints

BS 5194-3:1985

Specification for Dissecting Forceps

BS 5194-4:1985

Specification for Scissors, Shears and other Jointed Instruments

We operate our own UK production facilities equipped with ISO class 7 cleanrooms

DTR Medical

General & Vascular - Single-Use Instruments

DTR Medical Ltd

An Innovia Medical Company

17 Clarion Court, Enterprise Park, Swansea SA6 8RF UK

Telephone: +44 (0)1792797910

dtrinfo@innoviamedical.com

US 1.888.229.2875 | International +44 (0)1792797910

www.innoviamedical.com/single-use-instruments