

**SuperCut Stevens Tenotomy Scissors** Product Specification

**The Design Checklist**

- Ideal for delicate cutting and dissecting soft tissue
- One serrated blade and one sharp cutting blade allows a superior grip and cut
- Available with straight or curved blades
- Black handles make the instrument distinct



STSS2001



Code	Pack size	Description
DT- STSC2001	20 units	SuperCut Stevens Tenotomy Scissors Curved 115mm
DT- STSS2001	20 units	SuperCut Stevens Tenotomy Scissors Straight 115mm