

## Adson Forceps


## Product Specification

### The Design Checklist


- Narrow tips to effectively grip and hold delicate tissue in place
- Flexibility across a wide range of procedures



Adson Toothed Forceps  
- 2mm tip width



Adson Non-Toothed Forceps - 15mm tip length x 2mm tip width



Adson Forceps with Plateau- 1.4mm tip width



Fine Tip Adson Toothed Forceps - 22mm tip length x 1.5mm tip width



**DTR Medical**®

AN INNOVIA  
MEDICAL  
COMPANY

| Code         | Pack size | Description  |
|--------------|-----------|--|
| DT-ADS5001-S | 50 units  | Adson Non-Toothed Forceps 120mm, single packed                       |
| DT-ATS5001-S | 50 units  | Adson Toothed Forceps 120mm, single packed                           |
| DT-ATP2012-D | 20 units  | Adson Forceps 1x2 with plateau 120mm, double packed                  |
| DT-ATP2015-D | 20 units  | Adson Forceps 1x2 with plateau 150mm, double packed                  |
| DT-ATS2012F  | 20 units  | Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm, double packed |