


Beckmann Adenoid Curette

Product Specification

The Design Checklist

- Sharp and precise scoops to enable an effective Adenoidectomy
- Serrated handle to provide a secure grip whilst scraping Tissue



| Code | Pack size | Description |
|------------|-----------|---|
| DT-BAC2010 | 20 units | Beckmann Adenoid Curette, 10mm Curette Edge, 220mm, double packed |
| DT-BAC2013 | 20 units | Beckmann Adenoid Curette, 13mm Curette Edge, 220mm, double packed |
| DT-BAC2015 | 20 units | Beckmann Adenoid Curette, 15mm Curette Edge, 220mm, double packed |
| DT-BAC2017 | 20 units | Beckmann Adenoid Curette, 17mm Curette Edge, 220mm, double packed |