



Standard Bipolar Forceps

Product Specification

The Design Checklist

- Fine stainless steel tips with a range of sizes from 0.5 – 2mm
- Adjustable and universal 4mm plug that is tested to comply with the latest 'Flying Leads' UK legislation
- High quality cables and connectors deliver consistently repeatable results



Code	Pack size	Description
DT-SU770515044	20 units	Standard Bipolar Forceps 150 x 0.5mm, 3m, 4mm Plugs
DT-SU770520044	20 units	Standard Bipolar Forceps 200 x 0.5mm, 3m, 4mm Plugs
DT-SU771015044	20 units	Standard Bipolar Forceps 150 x 1mm, 3m, 4mm Plugs
DT-SU771520044	20 units	Standard Bipolar Forceps 200 x 1.5mm, 3m, 4mm Plugs
DT-SU772015044	20 units	Standard Bipolar Forceps 150 x 2mm, 3m, 4mm Plugs
DT-SU772020044	20 units	Standard Bipolar Forceps 200 x 2mm, 3m, 4mm Plugs