

**Hartmann Nasal Speculum**

**Product Specification**

**The Design Checklist**

- Hartmann Nasal Speculum featuring a strong stainless steel construction with smooth edges



Small - 32mm tip length -  
6.5mm tip width



Medium - 52mm tip length -  
6.5mm tip width



Large - 76mm tip length -  
6.5mm tip width



Code	Pack size	Description
DT-HNS2011-D	20 units	Hartmann Nasal Speculum, Small
DT-HNS2021-D	20 units	Hartmann Nasal Speculum, Medium
DT-HNS2031-D	20 units	Hartmann Nasal Speculum, Large