

Vannas Scissors

Product Specification

The Design Checklist

- Instrument of choice in Ophthalmic procedures such as Iridectomy and Capsulotomy
- Ideal for your most delicate practice, such as the removal of fine Sutures
- Variants are available in 85mm and 125mm
- Option of either Curved or Straight allowing flexibility
- Equipped with a ridged handle and a double leaf spring mechanism for self-opening action, ensuring enhanced control
- First-time sharpness which cut even at the tips
- Guaranteed sharp/sharp edges



Straight 85mm



Straight 125mm



Curved 85mm



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

Code	Pack size	Description
DT-VSS2081	20 units	Vannas Scissors Straight 85mm
DT-VSS2121	20 units	Vannas Scissors Straight 125mm
DT-VSC2081	20 units	Vannas Scissors Curved 85mm