

**Iris Scissors**

**Product Specification**

**The Design Checklist**

- Enhanced accuracy when it comes to delicate surgeries
- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips
- Available in both Straight and Curved variants
- Narrows to a fine and sharp point



Straight- 22mm tip length



Curved- 21mm tip length



Code	Pack size	Description
DT-ISC2110-D	20 units	Iris Scissors Curved 110mm
DT-ISC2111-S	20 units	Iris Scissors Curved 90mm
DT-ISS2110-D	20 units	Iris Scissors Straight 110mm
DT-ISS2111-S	20 units	Iris Scissors Straight 90mm