


Oral Rotating Biopsy Punch

Product Specification

The Design Checklist

- Designed specifically for use in oral procedures to obtain a biopsy from the front of mouth, gums, tongue, soft palate, and tonsil area.
- “Through cutting” jaw
- 360° rotation enables enhanced positioning and patient interaction
- Ergonomic handle designed for comfort
- Provides value-based benefits, improving procedure quality, patient safety, and overall cost savings
- Lightweight compared to alternatives
- UK manufactured



Code	Pack size	Description
DT-OBR1001	10 units	Oral Rotating Biopsy Punch 120mm
DT-OBRL1001	10 units	Oral Rotating Biopsy Punch Long 180mm