


Lacrimal Probes


Product Specification

The Design Checklist


- Double ended with 2 size options for adaptability during procedures
- Thin wire with rounded ends to effectively clear the passageway to the tear duct
- Blunt probe ends



0.7mm tip = 00



0.8mm tip = 0



0.9mm tip = 1



Code	Pack size	Description
DT-LCP2078	20 units	Lacrimal Probe, 0.7/0.8mm Diameter Wires
DT-LCP2089	20 units	Lacrimal Probe, 0.8/0.9mm Diameter Wires