

## Cusco Non-Insulated Specula without smoke tube

### The Design Checklist

- Strong metal construction with improved insulation performance for patient and user safety
- Smooth rounded design for improved patient comfort
- Robust self-retaining screw
- Small XXLong Metal Cusco Specula has been specially designed to improve the examination of post-menopausal women



*Cusco Specula Metal Small XXLong*



*Cusco Specula Metal Medium*

Code	Pack size	Description
DT-CMM200S	20 units	Cusco Specula Metal Small
DT-CMM200M	20 units	Cusco Specula Metal Medium
DT-CMM200ML	20 units	Cusco Specula Metal Medium Long
DT-CMM200L	20 units	Cusco Specula Metal Large
DT-CMM100SE	10 units	Cusco Specula Metal Small XXLong
DT-CMM100L	10 units	Cusco Specula Metal Large Long
DT-CMM100LE	10 units	Cusco Specula Metal Large Extra Long
DT-CMM100B	10 units	Cusco Specula Metal Broad

Width x Length (mm)	Insulated Smoke Tube	Metal Smoke Tube	Metal w/out Smoke Tube
<b>SMALL</b> 25 X 85mm	CSI200S	CSM200S	CMM200S
<b>MEDIUM</b> 30 X 85mm	CSI200M	CSM200M	CMM200M
<b>MEDIUM LONG</b> 30 X 95mm	CSI200ML	CSM200ML	CMM200ML
<b>LARGE</b> 35 X 95mm	CSI200L	CSM200L	CMM200L
<b>SMALL XXL</b> 25 x 135mm	-	-	CMM100SE
<b>LARGE LONG</b> 35 X 110mm	CSI100LL	CSM100LL	CMM100LL
<b>LARGE XL</b> 35 X 115mm	CSI100LE	CSM100LE	CMM100LE
<b>BROAD</b> 40 X 95mm	CSI100B	CSM100B	CMM100B