


Colposcopy LLETZ Ball Electrodes

Product Specification

The Design Checklist

- Single-use ensures no cross-infection risks or health and safety standards
- Gap between plastic and ball electrode for better access
- 152mm Ball Electrodes suitable for larger patient
- Non-Stick Teflon coating changes the heat dissipation and prevents tissue from sticking to the instrument



Code	Pack size	Description
DT- BET11303	10 units	Ball Electrode, 3mm diameter, Non-Stick Teflon Coated, 130mm length
DT- BET11305	10 units	Ball Electrode, 5mm diameter, Non-Stick Teflon Coated, 130mm length
DT- BET11503	10 units	Ball Electrode, 3mm diameter, Non-Stick Teflon Coated, 152mm length
DT- BET11505	10 units	Ball Electrode, 5mm diameter, Non-Stick Teflon Coated, 152mm length