

Colposcopy LLETZ Loop Electrodes

Product Specification

The Design Checklist

- Tungsten loops give ideal tension for greater control and are stiffer at smaller dimensions which prevents the wire from folding when cutting
- High cutting performance with reduced thermal damage to tissue
- Recognised colour coding with loop size printed on the handle to easily distinguish between the variants
- Universal shaft fits most electrosurgical pencils
- Single-use ensures first time sharpness and performance with no cross infection risks or health and safety standards



Code	Pack size	Description
DT- LEC1130505-PU	10 units	Loop Electrode, 5x5mm, 130mm Shaft (Purple)
DT- LEC1131005-LP	10 units	Loop Electrode, 10x5mm, 130mm Shaft (Light Purple)
DT- LEC1131010-YE	10 units	Loop Electrode, 10x10mm, 130mm Shaft (Yellow)
DT- LEC1131508-GY	10 units	Loop Electrode, 15x8mm, 130mm Shaft (Grey)
DT- LEC1131512-GN	10 units	Loop Electrode, 15x12mm, 130mm Shaft (Green)
DT- LEC1131515-OR	10 units	Loop Electrode, 15x15mm, 130mm Shaft (Orange)
DT- LEC1132010-RE	10 units	Loop Electrode, 20x10mm, 130mm Shaft (Red)
DT- LEC1132012-WH	10 units	Loop Electrode, 20x12mm, 130mm Shaft (White)
DT- LEC1132015-LB	10 units	Loop Electrode, 20x15mm, 130mm Shaft (Light Blue)
DT- LEC1132020-PI	10 units	Loop Electrode, 20x20mm, 130mm Shaft (Pink)
DT- LEC1132508-TA	10 units	Loop Electrode, 25x8mm, 130mm Shaft (Tan)