

# SAMARIT

safe patient moving & handling



Lightweight and comfortable, wherever you need it,  
in the clinic, home or away ...

## GLIDEBOARD

SEATED TRANSFER



**CURVED GLIDEBOARD**  
**THE GREAT**  
**THE BIGGEST**  
**SELFGlide**  
**TETRAglide**  
**ACCESSOIRES**



**SAMARIT**

Samarit transfer systems

Transfer from bed into the wheelchair or WC, wheelchair into the motorcar etc.



### THE BIGGEST

**Size:** 100 x 34 cm / 39 x 13 inch  
**Order No.** 470.0102

This board simplifies the transfer to the patient by the nursing personnel from the centre of the bed into the wheel-chair or WC/chair (even motorcar). It is manufactured of absolutely unbreakable material (HDPE), is acid resistant, easy to clean and to disinfect. A specially developed anti-slide adhesive (blue) ensures the positioning on the bed, wheelchair, etc.



### THE GREAT

**Size:** 80 x 33 cm / 32 x 13 inch  
**Order No.** 470.0080

This board simplifies the transfer from the edge of the bed into the wheelchair or WC, wheelchair into the motorcar etc. Its composition is the same as for "The Biggest". Sit-On Trousers (accessory) is used during the transfer without clothing, e.g. onto the WC-chair. It protects the skin of the patient and is fixed with Velcro. Two solid hand grips allow the nurse to carry out the transfer in a secure and controlled manner.



### CURVED GLIDEBOARD

**Size:** 68 x 23 cm / 27 x 9 inch  
**Order No.** 470.0012

This board serves for transfer bed-wheelchair-chair etc. We developed the shape together with therapists and several tens of thousands have been sold. The shape is Design Patent Protected. The material used corresponds to the other Glideboards (HDPE).



### SELFGlide

**Size:** 50 x 21 cm / 20 x 8 inch  
**Order No.** 470.0050

The Selfglide-Glideboard is primarily used by patients who can carry out the transfer themselves due to the strength in their arms. It can also be used in the transfer with the support of the Mini-care cushion by a nurse.



### TETRAGLIDE

**Size:** 63 x 23 cm / 25 x 9 inch  
**Order No.** 470.0863

The Tetraglide Glideboard was developed specially for tetraplegics but is also increasingly popular because of its new "Fishtail". It is characterised by the ease with which the patient can be sat upon, by the manner in which it can be inserted from both sides and is utilised for transfer bed-chair-WC, etc. The material corresponds to that of the other Glideboards.

Gap < 30 cm / 12 inch, SWL\* = 150 kg  
Gap < 20 cm / 8 inch, SWL\* = 180 kg



\*SWL = Safe Weight Load

### Accessories



### Sit-On Monofil Sliding-Towel

**Order No.** 470.0200

Sit-On Sliding-Towel permits easy sliding without clothes on to the Glideboards. Usually used by the patients themselves with the Curved, Selfglide and Tetraglide.  
Nylon for use to glide on Glideboards to bath, shower or toilet.



### Glide Belt

**Order No.** 470.2100 **Small**  
**Order No.** 470.2200 **Medium**  
**Order No.** 470.2300 **Large**

Helps to position the boards under heavy patients and protects the skin.

## GLIDEBOARD

Samarit transfer systems

- Unbreakable
- Sliding surfaces
- Non slide sticker
- Lightweight
- Left/right usable
- Instruction for use
- CE mark

**High quality standards**  
SAMARIT transfer aids have:

- Design patents
- 2 years warranty



[www.samarit.com](http://www.samarit.com)

Download your language from our website or ask your local distributor.



**Samarit Medical AG**

Gewerbestrasse 12

CH-8132 Egg, Schweiz

Tel.: +41 44 918 10 11 Fax: +41 44 918 24 39

e-mail: [info@samarit.com](mailto:info@samarit.com) - [www.samarit.com](http://www.samarit.com)